

We Think Better About Your Eyes!

EYE INSTRUMENTS CATALOGUE

MODREN RANGE OF EYE INSTRUMENTS



Manufacturer & Exporters of Surgical & Ophthalmic Instruments



Head Marala Road, Gohad Pur, Sialkot
- 51310, Pakistan.

Mob : 0092-305-6056248

E-Mail: info@dartersurgical.com

Web : www.dartersurgical.com

FORCEPS

Gill
Iris Forceps, curved
criss-cross serrated tips
1x2 teeth
FC-2001



Botvin
Iris Forceps
for peripheral iridectomy
serrated tips
1x2 teeth
FC-2002



Gills-Colibri
Forcep
very fine pointed tips
FC-2003



Colibri
Forcep, very delicate
1x2 teeth, 0.12mm
FC-2004



Harms-Colibri
Forcep, very delicate
1x2 teeth, 0.12mm
with tying platform
FC-2005



Dyson-Colibri
Forcep, very delicate
1x2 teeth, 0.12mm
extremely thin shanks
short tips with tying platform
FC-2006



FORCEPS

Colibri

Forceps, standard
1x2 teeth, 0.4mm
with tying platform
FC-2007



Polack

Corneal suturing forceps
1mm spread, very delicate
1x2 teeth, 0.12mm
FC-2008



Hofmann-polack

Corneal suturing forceps
with guide pin, 1mm spread
1x1 pierse teeth, 0.2mm

FC-2009



Pierse-colibri

Corneal forceps
very delicate
1x1 pierse teeth, 0.1mm

FC-2010



Pierse

Corneal forceps
very delicate
1x1 pierse teeth, 0.1mm
straight
curved

FC-2011



Castroviejo-Colibri

Forceps very delicate
1x2 teeth, 0.12mm

FC-2012



FORCEPS

Girards

Corneoscleral Forceps
1x2 teeth, 0.12mm
with tying platform
right handed
left handed
FC-2013



Maumenee-Colibri

Corneoscleral Forceps
very delicate
1x2 teeth, 0.12mm
FC-2014



Tennant-Colibri

Corneoscleral Forceps
very delicate
1x2 teeth, 0.12mm
FC-2015



Bores

Corneoscleral Fixation Forceps
U-shaped, 3mm spread
1x2 teeth at each tip,
0.12mm
FC-2016



Rubman Bores

Corneoscleral Fixation Forceps
U-shaped, 3mm spread
1x2 teeth at each tip,
0.12mm
FC-2017



Arrowsmith

Corneoscleral Fixation Forceps
with lock
U-shaped, 8mm spread
1x2 teeth at each tip,
0.12mm
FC-2018S
FC-2018A



FORCEPS

Dressing

Forceps, serrated
delicate
straight
curved

FC-2024S
FC-2014A



Tissue

Tissue Forceps
1x2 teeth
delicate
standard

FC-2025S
FC-2025A



Tissue

Forceps, 1x2 teeth
delicate
straight
curved

FC-2026S
FC-2026A



Lester

Fixation Forceps
straight
1x2teeth
2x3teeth

FC-2027S
FC-2027A



Barken

Fixation Forceps
1x2 teeth

FC-2028



FORCEPS

Katena
Corneal Fixation Forceps
with lock
U-shaped, 13mm spread
1x2 teeth at each tip,
0.12mm
straight
angled
FC-2019S
FC-2019A



Bishop-Harmon
Tissue Forceps
1x2 teeth
delicate
standard
FC-2020S
FC-2020T



Bishop-Harmon
Dressing Forceps
criss-cross serrated tips
FC-2021



Sauer
Fixation Forceps
1x2 angled teeth
FC-2022



Manhattan Eye & Ear
Fixation Forceps
1x2 angled teeth
with tying platform
FC-2023



FORCEPS

Farkas-Bracken
Fixation Forceps
1x2teeth, delicate
FC-2029



Kuhnt
Fixation Forceps
1x2 curved teeth
FC-2030



Elschnig
Fixation Forceps
1x2 angled teeth
FC-2031



Nugent
Superior Rectus Forceps
angled, 1x1 teeth
FC-2032



Troutman
Superior Rectus Forceps
angled, 1x2 teeth
FC-2033



Tennant-Troutman
Superior Rectus Forceps
curved shanks, 1x2 teeth
FC-2034



Graefe
Fixation Forceps
4.5mm wide jaws
without lock
with lock
FC-2035A
FC-2035B



FORCEPS

Tennant

Tying Forceps
with guide pin
extra delicate
smooth jaws with
6mm long platform
for 9-0 to 11-0 sutures
straight
curved



FC-2036S
FC-2036C

Harms

Tying Forceps
smooth jaws
for 7-0 to 9-0 sutures
straight
curved



FC-2037S
FC-2037C

Fine

Tying Forceps
smooth curved jaws
for 7-0 to 9-0 sutures
FC-2038



Barraquer

Cilia Forceps
lightly angled
smooth curved jaws
FC-2039



Katena

Cilia Forceps
slanted jaws
gold plated
FC-2040



Gradle

Cilia Forceps
smooth, squared off
positive gripping jaws
FC-2041



FORCEPS

Lambert
Chalazion Forceps
round, 8mm diameter
FC-2042



Hunt
Chalazion Forceps
round, 8mm diameter
FC-2043



Baird
Chalazion Forceps
oval, small
FC-2044



Francis
Chalazion Forceps
oval, medium
FC-2045



Desmarres
Chalazion Forceps
oval, large
FC-2046



Wies
Chalazion Forceps
grooved upper jaw
FC-2047



Lordan
Chalazion Forceps
triangular jaws
FC-2055



FORCEPS

**Jewelers
Forceps, delicate
FC-2048**



**Jewelers
Forceps, delicate
FC-2049**



**Jewelers
Forceps, extra delicate
FC-2050**



**Osher
Iris Forceps
very fine textured
jaws
FC-2051**



**Jewelers
Forceps, extra delicate angled
FC-2052**



**Jewelers
Forceps, curved
FC-2053**



**Suture Removal Forceps
Very fine pointed tips
with carbide impregnated jaws
FC-2054**



FORCEPS



Serrefine
Small, straight
FC-2055



Hartman
Mosquito Forceps
extra delicate
3 3/4" (95mm) long
straight
curved
FC-2056



Halsted
Mosquito Forceps
delicate
5 3/4" (145mm) long
straight
curved
FC-2057



Backhaus
Towel Clamp
3 1/2" (90mm) long
5 1/4" (135mm) long
FC-2058



Baby Jones
Towel Clamp
cross-action
2 1/4" (55mm) long
FC-2059



FORCEPS

Shepard II
Lens Holding Forceps
with nudge notch
for A/C IoLs
with lock
without lock
FC-2060



Jameson
Muscle Forceps
with slide lock
4teeth, child size
right
left
FC-2061R
FC-2061L



Jameson
Muscle Forceps
with slide lock
6teeth, adult size
right
left
FC-2062R
FC-2062L



Berke
Ptosis Forceps
with slide lock
longitudinally grooved jaws
20mm long jaws
27mm long jaws
FC-2063A
FC-2063B



Erhardt
Lid Forceps
FC-2064



SCISSORS



Castroviejo
Keratoplasty Scissors
angled to side
blunt tips
small blades
SC-2065



Castroviejo
Keratoplasty Scissors
angled to side
blunt tips
medium blades
SC-2066



McPherson-Westcott
Conjunctival Scissors
curved, blunt tips
small blades
SC-2067



Westcott
Tenotomy Scissors
curved, blunt tips
medium blades
SC-2068



Westcott
Tenotomy Scissors
curved, blunt tips
standards blades
right
left
SC-2069R
SC-2069L



EYE INSTRUMENTS CATALOGUE

Markers.Calipers

Anis

Radial Marker
40. mm center opening
12.5mm O.D.
3 blades 10 blades
4 blades 12 blades
6 blades 16 blades
8 blades

MC-2070A
MC-2070B



Ains

Suture Placement Marker
for corneal grafts, 8 radial blades
with 8mm center opening
MC-2071



Deitz

Incision Depth Gauge
for corneal incisions
0.5mm
0.55mm
MC-2072



Thornton

Limbal Incision Ruler
with marking lines
lines spaced 1mm apart*
MC-2073



Lindstrom

Small Incision Marker
for cataract incision widths of
3, 4, and 5mm at 1, 2, or 3mm
away from the limbus*
MC-2074



Castroviejo

Caliper, 20mm spread
straight
curved
MC-2075



Osher

Incision Caliper
for internal measurement of
cataract incisions from 1 to 10mm
MC-2076



SCISSORS



Castroviejo
Corneoscleral Scissors
small blades
right
left
MC-2077R
MC-2077L



Castroviejo
Corneoscleral Scissors
Medium blades
Right
Left
MC-2078R
MC-2078L



Troutman-Castroviejo
Corneoscleral Scissors
with stop, small blades
Right
Left
MC-2079R
MC-2079L



EYE INSTRUMENTS CATALOGUE

Vannas
Capsulotomy Scissors
sharp tips
angled to side
SC-2080



Gills-Vannas
Capsulotomy Scissors
extra thin long blades
sharp tips
medium blades
SC-2081



Gills-Vannas
Capsulotomy Scissors
extra thin long blades
sharp tips
SC-2082



Gills-Vannas
Capsulotomy Scissors
extra thin long blades
sharp tips
SC-2083



McPherson-Vannas
Iris Scissors
sharp tips
curved
SC-2084



Katena-Vannas
Scissors, extra thin
7mm long blades
sharp points
straight
SC-2085



FORCEPS



Shaepard
Lens holding Forceps
with serrated lowre jaw
gently curved shanks
without lock
with lock
FC-2086



Sinskey
Lens holding Forceps
delicate T-shaped lowre
jaw for 3-points fixation
FC-2087



Kraff
Lens holding Forceps
delicate curved jaws
to grasp edge of lens optic
FC-2088




Anis
Lens Inserting Forceps
very delicate curved tips
to grasp IOL haptics.
Angled shanks for
comfortable hand position
FC-2089



Clayman
Lens holding Forceps
delicate angled jaws
gently curved shanks
without lock
with lock
FC-2090



Phaco-Folder II
For folding and inserting
soft lols through a
small incision
FC-2091



EYE INSTRUMENTS CATALOGUE

Eye Scissors
Standard
pointed tips
straight
curved
SC-2092S
SC-2092C



Eye Scissors
Standard, larg rings
pointed tips
straight
curved
SC-2093S
SC-2093C



Stitch Scissors
Ribbon type
needle points
light curved
SC-2094



Spatulas



Barraquer
Corneoscleral brush
Sable with plastic handle
SP-2095



McPherson
Iris Spatula
1mm wide
malleable stainless steel
SP-2096



Culler
Iris spatula
Sterling silver
Blade
1mm wide
2mm wide
SP-2097A
SP-2097B



Wheeler
Double-Ended Spatula
Lightly curved
SP-2098



Barraquer
Iris Spatula
Round, smooth
0.5mm dia.
0.25 mm dia.
SP-2099



Gimbel
Nucleus Spatula
Angled, round, smooth
0.25mm dia.
SP-2100



Knolle
Lens Nucleus Spatula
Malleable stainless steel
Flat, 0.5mm wide
SP-2101



Nucleus Manipulators

Bechert
Lens Nucleus Rotator
Blunt forked tip
Angled at 7mm
NM-2102



Jaffe-Bechert
Lens nucleus Rotator
Blunt forked tip
Angled at 9mm
NM-2103



Clayman
Nucleus Rotator
0.5mm wide with
V-notched tip
NM-2104



Drysdale
Nucleus Manipulator
Paddle shaped tip
NM-2105



Pisacano
Nucleus Rotator
Flat blunt tip
NM-2106



Siepser
Endocapsular Controller
Trumpet shaped tip
NM-2107



Tennant
Nucleus Manipulator
Ball shaped tip
NM-2108



Spatulas.loops



Elsching
Cyclodialysis Spatula
Flat, 1mm wide blade
SL-2109



Castroviejo
Synechiaie Spatula
semi-sharp tips
0.5mm wide
double-end
SL-2110



Castroviejo
Cyclodialysis Spatula
double-end
0.25mm wide
0.5mm wide
SL-2111A
SL-2111B



Kimura
Platinum Spatula
SL-2112



Lewis
Lens Loop
Small
Large
SL-2113S
SL-2113L



New Orleans
Lens Loop
SL-2114



Wilder
Lens Loop
serrated
SL-2115



Punches.polishers

Gass
Corneoscleral punch
1.5mm dia. Tip
rotatable for
Convenient positioning
PP-2116



Kelly
Descemet's
Membrane Punch
1.0mm dia.tip
PP-2117



Meyerhoefer
Chalazion Curette
Size 0, 1.5mm dia.
Size 1, 1.75mm dia.
Size 2, 2.25mm dia.
Size 3, 2.5mm dia.
Size 4, 3.5mm dia.



PP-2118A PP-2118B PP-2118C PP-2118D

Gills
Capsule Polishing Curette
1.75mm dia. Cup
with sandblasted base
gentel backward angle
PP-2119



Nightingale
Capsule Polishing Curette
1.75mm dia. Open ring
with semi-sharp edge
PP-2120



Anis
Capsule Scrubber
Gently curved flexible shank
With roughened tip
Ball tip
Disc tip
Reverse ball tip



PP-2121A PP-2121B PP-2121C

Spatulas.manipulators



Lippman
Sideport Nucleus Spatula
0.5mm dia. Round
blunt forked tip
SM-2122



Girard
Synechiae Spatula
Semi-sharp tip
0.5mm wide
SM-2123



Obstbaum
Lens Spatula
Angled blade
0.25mm wide
SM-2124



Hirschman
Lens Spatula
Flat blade, 0.5mm wide
With small notches for
Manipulating IOLs
SM-2125



Jaffe
Lens Spatula
Flat, smooth
0.5mm wide
SM-2126



Lester
IOL Manipulator
For Lenses with or
Without hole
Straight, with guard
angled
SM-2127S
SM-2127A



EYE INSTRUMENTS CATALOGUE

Loops

Kirby
Lens Loop
And Muscle Hook
L -2128



Gills-Welsh
Lens Loop and Probe
with serrations
L -2129



Morrison
Lens Loop
And Probe
Very thin
L -2130



Simcoe
Lens Nucleus Loop
Fine serrations, strong curve
Double-end
L -2131



Kansas
Lens Loop, 3mm wide
Serrated, for small incision
Nucleus fragment
Removal
L -2132



Keener-Arit
Lens Loop
3mm wide, serrated
strong curve for
small incision nucleus
fragment removal
L -2133



Keener
Nucleus Divider
With two pre-formed
Snare wire loops
Complete
Wire loops(pk/3)
L -2134



Spoons.irisretractors



Bunge
Evisceration Spoon
Small
Large
SI-2134S
SI-2134L



Wells
Enucleation
Spoon
SI-2135



Rizzuti
Iris Retractor
And Lens
Expressor
2mm wide lip
SI-2136



Azar
Iris Retractor
Autoclavable plastic
Left
Right
SI-2137L
SI-2137R



Graether
Collar Button
Micro iris retractor
Straight, with guard
angled
SI-2138S
SI-2138A



Nucleus dividers.spoons

Kansas
Nucleus Trisector
For manual small incision surgery
For right handed
Surgeons
For left handed
Surgeons
NS-2139



Kansas
Nucleus Bisector
For manual small incision surgery
NS-2140



Kansas
Nucleus Vectis, solid
Used as a "cutting block"
For Kansas nucleus
Bisector and trisector
For right handed surgeons
For left handed surgeons
NS-2141



Culler
Lens Spoon
Sterling Silver
NS-2142



Barraquer
Lenticular Spoon
For corneal button transfer
NS-2143



Paton
Spatula and Spoon
For corneal button transfer
NS-2144



Iris hooks



Graether
Collar Button
Micro iris retractor
And posterior Capsule polisher
With side irrigation
Straight
Angled
IH-2145S
IH-2145A



Bonn
Micro Iris Hook
Blunt, round handle
With protective Guard
IH-2146



Bonn
Micro Iris Hook
Blunt, flat handle
IH-2147



Jaffe-Knolle
Iris Hook, very blunt
Small, angled
IH-2148



Knolle
Iris Hook, very blunt
Small, straight
IH-2149



Shepard
Iris Hook
Small, blunt
IH-2150



IOL hook

Sinskey
Lens Manipulating Hook
With round handle
0.25mm dia. Blunt tip
straight, with guard
angled
IH-2151S
IH-2151A



Sinskey II
Lens manipulating Hook
With round handle
0.2mm dia. Blunt tip
straight
angled
IH-2152S
IH-2152A



Sinskey II
Lens Manipulating Hook
With round handle
0.2mm dia. Blunt tip
straight, with guard
angled
IH-2153S
IH-2153A



Lewcky
Lens manipulating Hook
0.15mm dia. Blunt tip
angled
IH-2154



Hirschman
Iris Hook
Blunt, angled
IH-2155



Hooks



Maltzman-Fenzi
 Lens Manipulating Hook
 0.15mm dia. V-shaped tip
 round handle
 straight, with guard
 angled
H -2156



Galand
 Lens Manipulating Hook
 angled
H -2157



Jaffe-Maltzman
 Lens Manipulating Hook
 0.15mm dia. V-shaped tip
 angled
H -2158



Osher
 Y-Hook
 for manipulating IOLs
 straight, with guard
H -2159



Fixation Forceps
 double, sharp
 small
 large
H-2160S
H-2160L



Twist Hook
 for scleral fixation
 right
 left
H-2161R
H-2161L



Anis
 Scleral Twist Fixation Cone
 2.5mm dia. With three
 teeth, with guard
H -2162



EYE INSTRUMENTS CATALOGUE

Hooks

Kuglen
Iris Hook
And Lens Manipulator
Push-pull model
Straight
Angled
H-2163S
H-2163A



Kuglen
Iris Hook
And Lens Manipulator
Clover leaf model
Straight
H -2164



Kuglen
Iris Hook
And Lens manipulator
Push-pull model
With round handle
Straight, with guard
Angled
H-2165S
H-2165A



Rentsch "Boat Hook"
Lens Manipulator
This multi-purpose
Instrument may be used as
An iris/capsule retractor or
Pusher/manipulator for IOLs
With and without holes
Straight
With guard
H -2166



Fenzl
Lens Manipulating Hook
0.15mm dia. T-shaped tip
straight,
with guard
angled
H-2167S
H-2167A



Needle Holders



McPherson
Needle Holder
very delicate jaws
straight, with lock
straight, without lock
curved, with lock
curved, without lock
NH-2168A NH-2168C
NH-2168B NH-2168D



Jaffe
Needle Holder
delicate conical jaws
straight, with lock
curved, with lock
curved, without lock
NH-2169A
NH-2169B
NH-2169C



Castroviejo
Needle Holder
delicate jaws
straight, with lock
straight, without lock
curved, with lock
curved, without lock
NH-2170A
NH-2170B
NH-2170C
NH-2170D



Needle Holders

Castroviejo
Needle Holder
standard jaws, with lock
straight
curved
NH-2171A
NH-2171B
NH-2171C



EYE INSTRUMENTS CATALOGUE

Kalt
Needle Holder
straight smooth jaws
delicate
NH-2172



Drews
Needle Holder
extra delicate jaws
straight, with lock
straight, without lock
curved, with lock
curved, without lock
NH-2173A
NH-2173B
NH-2173C
NH-2173D



Barraquer
Needle Holder
short model
very delicate jaws
straight, with lock
straight, without lock
curved, with lock
curved, without lock
NH-2174A NH-2174C
NH-2174B NH-2174D



Needle Holders



Anis

Needle Holder
extra delicate jaws
straight, with lock
straight, without lock
curved, with lock
curved, without lock
NH-2175S
NH-2175C



Katena

Needle Holder
extra delicate jaws
straight, with lock
straight, without lock
curved, with lock
curved, without lock
NH-2176A
NH-2176B
NH-2176C
NH-2176D



McInyre

Fish Hook Needle Holder
extra delicate jaws
straight, with lock
NH-2177



Needle Holders

Troutman
Needle Holder
delicate, strongly curved jaws,
with lock
without lock
NH-2178A
NH-2178B



Barraquer
Needle Holder
very delicate curved jaws
with lock
without lock
NH-2179A
NH-2179B



Barraquer
Needle Holder
delicate jaws
straight, with lock
straight, without lock
curved, with lock
curved, without lock
NH-2180A
NH-2180B
NH-2180C
NH-2180D



Barraquer
Needle Holder
short model
standard jaws, curved
with lock
without lock
NH-2181



Speculums



Barraquer
Wire Speculum
 small
 large
 large, heavy wire
 SP-2182S
 SP-2182L



Barraquer
Wire Speculum
 with solid blades
 small
 large
 SP-2183S
 SP-2183L




Grandon-Barraquer
Wire Speculum
 with suture notches or holes
 to obtain additional traction
 for the upper lid
 wire blades
 solid blades
 SP-2184W
 SP-2184S



Feaster-Barraquer
Wire Speculum
 wide, solid blades
 SP-2185



McIntyer
Wire Speculum with
 gently curved spring.
 V-shaped blades help to retain
 disposable drape under lids
 close wire
 open wire
 SP-2186




EYE INSTRUMENTS CATALOGUE

Speculums

Kartz-Barraquer
Wire Speculum
with open blades
for unobstruted access
to the cataract incision
with phaco or I/A tip
heavy wire, large
standard wire, large
standard wire, small
SP-2187L
SP-2187S



Kartz-Barraquer
Wire Speculum with wings
for retracting surgical
drape while allowing
unobstruted access
to the cataract incision
with phaco or I/A tip
standard wire, small
SP-2188W
SP-2188S



EYE INSTRUMENTS CATALOGUE

Tennant
Lid Retractor
may be clamped to surgical
drape with suture or elastic band
pair
SP-2189



Sauer
Lid Speculum
for infants
for premature infants
SP-2190

Infant
wire Speculum
with gently curved springs
4mm wide blades
SP-2191



Cook
Eye Speculum with looking screw
solid,small #1
solid,medium #2
solid,large #3
fenestrated, medium #2
fenestrated, large #3
SP-2192S
SP-2192M SP-2192L

Speculums



Castroviejo
Eye Speculum
with finely controlled opening
and closing mechanism
small #1
SP-2193



Castroviejo
Eye Speculum
with finely controlled opening
and closing mechanism
medium #2
SP-2194



Fishkind-castroviejo
Lid Speculum with open blades
for use with phaco I/A tip
SP-2195



McPherson
Eye Speculum
solid blades with
vertical suture posts
SP-2196



Speculums

Williams
Eye Speculum
open blade, small
SP-2197



Williams
Eye Speculum
open blade, large
SP-2198



Lancaster
Eye Speculum
solid blade
SP-2199



Lester-Burcs
Eye Speculum
solid blade with
flanges for retaining
eye lashes
SP-2200



Speculums



Meaumennee-park

Eye Speculum
with canthus hook,
blade locking screws may
be positioned vertically
or horizontally
solid blades
fenestrated blades
SP-2201



Kuglen

Eye Speculum
with canthus hook,
solid blades with
swivel ball joints
and locking screws
SP-2202



Guyton-park


Eye Speculum
with vertical suture posts
and canthus hook
solid blades
fenestraeed blades
SP-2203




Smirmaul
Eye Speculum
individually adjustable
blade for elevating lids
solid blades
wire blades
SP-2204



Schott
Eye Speculum
with rotatable adjustable
for precise control of
blades elevatin
wire blades
SP-2205



Flieringe
Scleral Fixiation Ring
set of 8 sizes(15-22mm)
14mm dia.
15mm dia.
16mm dia.
17mm dia.
18mm dia.
19mm dia.
20mm dia.
21mm dia.
22mm dia
SP-2206



Retractors



Jaeger
Lid plate
stainless steel
RE-2207



Wright
Fascia Needle
with 1mm x 6mm
oval hole
14cm long
RE-2208



Knapp
Lacrimal Retractor
4-prong
blunt
sharp
RE-2209



Stevenson
Lacrimal Sac Retractor
3 x 3 blunt prongs
self-retaining
adjustable spread
RE-2210



Stevenson
Lacrimal Sac Retractor
solid blades with
serrated edge
self-retaining
adjustable spread
RE-2211



Agricola
Lacrimal Sac Retractor
3 x 3 sharp prongs
RE-2212



Retractors

Schepens
Froked Orbital Retractor
SP-2213



Arruga
Globe Retractor
SP-2214



Jaffe
Wire Lid Retractor
may be clamped to
surgical drape with
suture or elastic band
pair,large
pair,small
SP-2215



Desmarres
Lid Retractor
solid blade
size 0
size 1
size 2
size 3
SP-2216A
SP-2216B
SP-2216C
SP-2216D



Jaeger
Lid Plate
clear plastic
SP-2217



Knives



Zieglar
Knife Needle
3mm
5mm
KN-2218



Girard-swan
Knife Needle
1.5mm cutting edge
5mm
KN-2219



Wheeler
Discission Knife
18mm cutting edge
KN-2220



Emery
cystotome
double-edge
for push and pull cutting
blade, straight
blade, angled
chuck handle
KN-2221



Stiletto Blade
For micro-vitro-retinal
stab incisions, lancet
shaped, 19 gauge
1.2mm wide
blade
chuck handl
KN-2222



Paufique
Graft Knife
angled tip
curved cutting edge
KN-2223



Knives.spuds

Tooke
Coeneal Knife
straight blade with
curved cutting edge
KS-2224



Gill
Coeneal Knife
curved cutting edge
straight
curved
KS-2225



Golf-Club
Foreign Body Spud
angled, lightly rounded tip
KS-2226



Francis
Foreign Body Spud
angled, flag shaped tip
KS-2227



Ellis
Foreign Body Spud
very thin, lightly
curved and rounded tip
KS-2228




Spud & Gouge
In reversible handle
for removing foreign bodies
KS-2229



Wills Eye Hospital
Foreign body spud & bar
delicate curved spud on
one end and small barr
for rust ring reoval
on opposite end
KS-2230



Cannulas



West
Irrigating Cannula
Blunt tip, 23 gauge
Side Opening
KS-2231



Kelman
Cyclodialysis cannula
spatulated tip
side opening, 19 gauge
left opening
Right opening
KS-2232



Air injection
Cannula, angled 5mm from tip
30 gauge
27 gauge
KS-2233



Sheets
Cyclodialysis Cannula
spatulated tip
one front and two side
openings, 19 gauge
KS-2234



Viscoelastic
Aspirating cannula
Angled 10mm from tip
22 gauge
19 gauge
KS-2235



Troutman
Alpha-Chymotrypsin Cannula
delicate olive tip
25 gauge
KS-2236



Rainin
Air injection cannula
Spatulated tip
Angled at 7mm
30 gauge
27 gauge
KS-2237



Randolph
Cyclodialysis cannula
Curved spatulated tip
Angled at 13mm
21 gauge
19 gauge
KS-2238



Wfelsh
Flat Olive tip Cannula
25 gauge
straight
left
right
KS-2239S
KS-2239L
KS-2239R

Blade breakers

Aldgerbrush
Corneal Rust
Ring Remover
battery powered
complete
extra burr 0.5mm
extra burr 1.0mm
KS-2240



Castroviejo
Blade Breaker
and holder
with sliding lock
KS-2241



Castroviejo
Blade Breaker and holder
short model
pointed jaws
squared jaws
KS-2242



Swiss Model
Blade Breaker
and holder
spherical tip
KS-2243




Katena
Double Edge Blades
for blade breakers
pk/10 non-sterile
pk/50 sterile
KS-2244





Elliot
Scleral Trephine
1.0mm dia.
1.5mm dia.
2.0mm dia.
Handel
KS-2245



I-A Cannulas




Simcoe
 Irrigating-Aspirating Cannula
 Original model, with aspiration
 Through tubing hub and irrigation
 Through luer-lock hub.23 gauge
 Thin wall, 15mm long curved
 Tubes, with 0.3mm dia. Aspirating
 Port facing into curve and slightly
 Shorter irrigating tip with 45°
 Angled front opening. Designed
 For true Simcoe I/A technique.
IA-2246



Swiss
 Irrigating-Aspirating cannula
 for aspiration through tubing hub
 and irrigation through luer-lock
 hub.23 gauge thin wall, 15mm long
 curved tubes,0.3mm dia. Port on left
 side of aspirating tube. Slightly
 angled front opening.
IA-2247



Reversed simcoe
 Irrigating-Aspiratin Cannula
 Same as K7-4300 but tubes reversed.
 For surgeons who like the simcoe I/A
 Cannula, but prefer to use a more
 conventional technique of aspirating
 through the luer-lock hub cannula .
IA-2248



Ambrosio
 Irrigating-Aspirating
 Cannula For aspirating
 through tubing hub and
 Irrigation through luer-lock
 hub.23 gauge thin wall,
 15mm long, curved,
 front openings.
 smooth tip
 rough tip for polishing
IA-2249


EYE INSTRUMENTS CATALOGUE

Cystotomes

EYE INSTRUMENTS CATALOGUE




Kelman
Cystotome
With sharp point and
Sharp cutting edge
22 gauge
IA-2250



Knolle-Kelman
Cystotome
small blade
with sharp point and
sharp cutting edge
22 gauge
IA-2251



Kratz
Cystotome
With sharp point and
Sharp cutting edge
Vaulted, 22 gauge
IA-2252



McIntrye
Cystotome
angled 10mm from tip
with sharp point, 23 gauge
IA-2253



Kelman
Cystotome
With sharp point and
Blunt blade, 22 gauge
IA-2254



Lieppman
Micro Cystotome
extra small, with sharp point
and sharp cutting edge
22 gauge
IA-2255



Kelman
Double Blade Cystotome
With sharp points and
Blunt blades, 22 gauge
IA-2256



Graether
Collar Button
micro iris retractir
and posterior capsule polisher
with side irrigation
straight
angled
IA-2257



Double-Edge
Cystotome
With sharp point and
sharp cutting edges
22 gauge
IA-2258



Shammas
Posterior Capsule
Polishing Curette
1mm diameter cup with
sandblasted base
22 gauge
IA-2259

K-Blades. handels



BH-2274A Miniature Blade
Sterile, disposable
#63 (pk/6)
#64 (pk/12)
#65 (pk/12)
#67 (pk/12)
#69 (pk/6)
#69B (pk/6)



Lamellar Blade
sterile, disposable (pk/3)
#66-60
#66F
BH-2273A
BH-2273B



Super Blade
sterile, disposable (pk/12)
#76-15,15*
#76-30,30*
#76-45,45*
BH-2275A
BH-2275B
BH-2275C

Chuck Handle
#3H hexagon 10cm
#3K round 10cm
#3KA round 7.5cm
#3KB round 5.0cm

BH-2276A
BH-2276B
BH-2276C
BH-2276D



Stiletto Knife
For Micro-Vitreo-Retinal
Stab incisions
Sterile, disposable (pk/3)
20 gauge
19 gauge
BH-2277



Bard-Parker Handle
For disposable slot blades
#3
#9
BH-2278A
BH-2278B



Bard-Parker Blade
Sterile, disposable (pk/50)
#11
#15
BH-2279A
BH-2279B




Capsule Polishers




Drews
Posterior Capsule Polisher
Lightly angled, 2mm I.D loop
With sandblasted base
0.2mm irrigating port facing
into circle, 25 gauge
IA-2263



Jensen
Posterior Capsule Polisher
carbide impregnated on
Posterior surface of olive tip
front opening, gently curved
25 gauge
23 gauge
IA-2264



Kraff
Posterior Capsule
Polishing Curette
1.75mm diameter cup
with semi-sharp edge.
Irrigating port in center
Of cup, 25 gauge
IA-2265



Jensen
Capsule Polisher
carbide impregnated
All around olive tip
front opening, gently curved
25 gauge
23 gauge
IA-2266



Kratz
Scratcher
Front opening with
Carbide impregnated tip
Lightly angled, 22 gauge,
IA-2267



Knolle
Posterior Capsule Polisher
2mm angled tip, front opening
finely sandblasted
30 gauge
23 gauge
IA-2268




Rubman-katzin
Posterior Capsule Polisher
Carbide impregnated tip
With smooth closed end.
Irrigating side port facing
Into curve, 22 gauge
IA-2269



Simcoe
Cortex Extractor and Capsule Polisher
side opening facing into curve
23 gauge, sandblasted
13mm
21mm
IA-2270



Shapiro
Posterior Capsule Polisher
Carbide impregnated tip
With dual 0.3mm sideports
curved. 23 gauge
IA-2271




Hydro-Dissection
Cannula, 25 gauge
Angled flat tip
IA-2272


Disposable cannulas




Bishope-Harmon
Irrigating Cannula
Flat tip, angled at
8mm from the tip
sterile, disposable
19 gauge
DC-2280



Simcoe
Irrigating-Aspiration
Cannula, 23 gauge thin wall
0.3mm aspiration port
with tubing and adaptor
sterile, disposable
DC-2282




Hydrodissection
Cannula, flat tip
Angled at 7mm
From the tip
Sterile, disposable
25 gauge
DC-2281



Kartz
Capsule Polisher
end opening, angled at
7mm from the tip
sterile, disposable
23 gauge
21 gauge
DC-2284



Hydrodissection
cannula, U-Shaped
with flat tip
sterile, disposable
25 gauge
DC-2283



Kartz
Capsule Polisher
close tip with side port, angled
at 7mm from the tip
sterile, disposable
23 gauge
21 gauge
DC-2286



J-Shaped Cannula
Sterile, disposable
Straight, 25 gauge
angled, 25 gauge
straight, 23 gauge
angled, 26 gauge
DC-2285



Simcoe
Cortex Extractor
23 gauge X 13mm
0.3mm aspiration port
sterile, disposable
DC-2287



Silicone
Tip Cannula
angled at 7mm from
the tip, 27 gauge
sterile, disposable
DC-2288

Disposable instruments

Bipolar Forceps
Jeweler's Type
pointed tips, straight
plain
insulated
19 gauge
DI-2289



Bipolar Forceps
McPherson Type
lightly blunted tips
straight
DI-2290



Bipolar Forceps
McPherson Type
lightly blunted tips
angled
DI-2291




EYE INSTRUMENTS CATALOGUE

Bipolar
Coaptation Forceps
pointed tips
DI-2292



Bipolar Cord
disposable(10/box)
reusable
DI-2293



Cannulas



Mcintyre-Binkhorst
Irrigating Cannula
smooth blunt tip
front opening, 26 gauge
right
straight
left
C-2260R
C-2260S
C-2260L



Gimbel
Hydro-dissection Cannula
smooth blunt tip
front opening, 30 gauge
C -2263



Knolle-Pearce
Irrigating Vectis
dual fluid flow through
three front opening ports
23 gauge loop
C -2264



Drows
Irrigating Vectis
dual fluid flow through
one front and two backward
opening ports, 23 gauge loop
C -2266



Sheets
Irrigating vectis
dual fluid flow through
three large front opening
ports, serrated 20 gauge loop
C -2261



Hofmann
Nucleus loops
3mm wide, serrated
extra delicate
non-irrigating
C -2262



Simcoe
Irrigating Lens Nucleus Loop
two front opening ports
finely serrated, strong curve
25 gauge
right
left
C-2265R
C-2265L



Anis
Irrigating Vectis
dual fluid flow through
two front opening ports
23 gauges heart shaped loop
enables surgeon to retract
iris and capsular flap
while irrigating
C -2267

cannulas



Nishi
Nucleus Delivery Cannula
flat insertion plate with
central irrigating port.
Restraining fin limits
depth of insertion and
may be used to apply gentle
pressure to the sclera.
Tip of insertion plate is
placed under the nucleus,
while continuous irrigation
helps to float the nucleus
forward across the plate,
reducing zonular stress.
Right
left
C -2268



Saiz
Irrigating Cassule Retractor
20 gauge, straight
C - 2269



Kratz
Lens Manipulating Cannula
0.2mm diameter tip
front opening, 30 gauge
C -2270




Thurmond
Nucleus Irrigating Cannula
and Iris Retractor, front opening
smooth downward curved tip
19 gauge
C -2271



Osher
Lens Manipulating Cannula
with sideport and pigtail tip
24 gauge
C -2272




Nishi
Irrigating Capsule Rectors
for hydro-dissection of the
lens nucleus
C -2273



Kuglen
Irrigating Lens Manipulator
disc shaped tip with 0.25mm
diameter peg protruding
from its center, 25 gauge
straight
angled
C-2274S
C-2274A

Cataract set



PLANNED EXTRA-CAPSULARA CATARACT EXTRACTION

CS-2275

- | | |
|---|---|
| <ul style="list-style-type: none"> BSS Irrigator Towel Clamp Mosquite Forceps Serrefine (2) Wire Speculum Super Sharp Blade Blade Handle Corneal Section Scissors, right Corneal Section Scissors, Left Westcott Scissors, blunt Stitch Scissors Capsulotomy Scissors Irais Scissors Caliper Conjunctive Forceps Superior Rectus Forceps Colibri Forceps 0.12mm Diamond Step Knife Intraoperative Keratometer Tissue Forceps Suturing Forceps, 0.12mm Tying Forceps, Stright Tying Forceps, angled Kelman-mcPherson Forceps | <ul style="list-style-type: none"> Utrata Capsulorhexis Forceps IoL Forceps IoL Manipulator Kuglen Iris Hook Cyclodialysis Spatula Needle Holder, micro Needle Holder, strright Lens Loops & Probe Thornton Limbal Ruler Kansas Nucleus Trisector Kansas Nucleus Bisector Kansas Vectis, solid Kansas Nucleus Loop Kansas Nucleus Divider Loop for Nucleus Divider Kansas Nucleus Fragment Forceps Viscoelastic Aspirating Cannula, 22ga. Air Injection Cannula, 30ga. Viscoelastic Injection Cannula, 26ga. Cystotome Capsule Polisher I/A Cannula Sterilizing Case(not illustrated) K-Spong Spears(not illustrated) |
|---|---|


Enucleation set



ENUCLEATION

Speculum
Tissue Forceps
Mosquito Forceps
Tenotomy Scissors
Enucleation Scissors
Muscle Hook
Sterilizing Case(not illustrated)
ES-2276

Radial keratotomy set




RADIAL KERATOTOMY

RKS-2276

Wire Speculum, solid
Diamond Micrometer, double edge
Blade Gauge
Marking Pen
Thornton Fixation Ring
Fixation Forceps, straight
Fixation Forceps, angled
Visual Axis Marker
Optic Zone Maker, 3.0mm
Optic Zone Maker, 3.25mm
Optic Zone Maker, 3.5mm
Optic Zone Maker, 3.75mm

Optic Zone Maker, 4.0mm
Optic Zone Maker, 4.25mm
Optic Zone Maker, 4.5mm
Optic Zone Maker, 5.0mm
Radial Marker 6 Line
Radial Marker 8 Line
Incision Spreading Forceps
Incision Depth Gauge
Irrigating Cannula, 27ga
Sterilizing Case (not illustrated)
K-Spong Spears (not illustrated)

Lid set




LID SURGERY

LS-2277

- BSS Irrigator
- Towel Clamp(2)
- Mosquito Forceps, curved(2)
- Tenotomy Scissors
- Utility Scissors
- Miniature Blade#69
- Blade Handle
- Westcott Scissors, blunt
- Stitch Scissors
- Ptosis Forceps
- Needle Holder, delicate
- Kalt Needle Holde
- Lid Plate
- Caliper
- Lid Retractor #1
- Lid Retractor #2
- Lid Retractor #3

- Muscle Hook#1
- Muscle Hook#1
- Fixation Hook, sharp
- Chalazion Curette #1
- Chalazion Curette #2
- Chalazion Curette #3
- Chalazion Curette #4
- Fixation Frceps
- Conjunctive Forceps
- Tissue Forceps
- Fixation Forceps,angled
- Suturing Forceps
- Tying Forceps, curved
- Chalazion Forceps, round
- Chalazion Forceps,ovel
- Lid Forceps
- Sterilizing Case(not illustrated)
- K-Spong Spears(not illustrated)

Lacrimal set




LACRIMAL SURGERY

LS-2278 (Dacryocystorhinostomy & Intubation)

- BSS Irrigator
- Towel Clamp(2)
- Mosquito Forceps, straight(2)
- Mosquito Forceps, curved(2)
- Kerrison Rongeur
- Nasal Speculum
- Bone Rongeur
- Mallet
- Maniature Blade #65
- Blade Handle
- Self Retaining Retractor
- Lacrimal Cannula, straight
- Lacrimal Cannula, curved
- Lacrimal Dilator
- Periosteal Elevator
- Lacrimal Chisel, straight
- Lacrimal Chisel, straight
- Pigtail Probe
- Retrieving Hook
- Intubation set
- Lacrimal probe set
- Dressing Forceps
- Tissue Forceps
- Utility Forceps
- suturing Forceps
- Rake Retractor(2)
- Stevens Scissors
- Utility Scissors
- Needle Holder
- Stitch Scissors
- K-Spong Spears(not illustrated)

Enucleation set



ASTIGMATIC KERATOTOMY

ES-2279

Wire Speculum, solid
 Diamond Micrometer, double edge
 Blade Gauge
 Marking Pen
 Thornton Fixation Ring
 Fixation Forceps, double
 Fixation Forceps, 0.12mm
 Visual Axis Marker
 Degree Gauge

Axis Marker
 T-Incision Marker
 Optic Zone Maker
 Thornton Corneal Ruler
 Arcuate Incision Marker
 Irrigating Cannula, 27 ga.
 Incision Spreading Forceps
 Incision Depth Gauge
 Irrigating Cannula, 27ga
 Sterilizing Case (not illustrated)
 K-Spong Spears (not illustrated)

Sterilizing Case




Silicone Sterilizing Case
Venting Holes
26cm x 15cm x 2.5cm
SC-2280



Silicone Sterilizing Case
Venting Holes
26cm x 24cm x 2.5cm
SC-2281




Silicone Sterilizing Case
Venting Holes
Full silicone finger mat top & bottom
26cm x 24cm x 2.5cm
accomadetes 16 instruments
SC-2282




Silicone Sterilizing Case
Venting Holes
26cm x 15cm x 2.5cm
accomadetes 10 instruments
SC-2283

EYE INSTRUMENTS CATALOGUE


Titanium Tissue Forceps




**ADSON
BD-30T**
120 mm




**ADSON
BD-31T**
120 mm




BD-141T
120 mm
9x9 Teeth




BD-09T
130 mm




BD-75T
145 mm



BD-16T
160 mm





**COLLEGE
TT-01**
150 mm





**MERIAM
TT-02**
160 mm

Titanium Tissue Forceps




1/1

GERALD
MF-101T
180 mm





1/1

GERALD
MF-102T
180 mm



1/1


GERALD
MF-103T
180 mm



1/1

GERALD
MF-104T
180 mm

Titanium Tissue Forceps




BD-32T/BD-33T
100 mm

BD-35T/BD-36T
100 mm

BD-60T/BD-62T
110 mm

BD-42T/BD-43T
110 mm



BD-285T
90 mm

BD-270T
130 mm


BD-341T
115 mm

**SEMKEN
BD-46T**
130 mm

Titanium Tissue Forceps



Mikro-Pinzette
Ø 8 mm



0,4 mm
1/1
MF-130T
150 mm

0,4 mm
1/1
MF-131T
150 mm

0,8 mm
1/1
MF-132T
150 mm

0,8 mm
1/1
MF-133T
150 mm

0,4 mm
1/1
MF-134T
180 mm

0,4 mm
1/1
MF-135T
180 mm

0,8 mm
1/1
MF-136T
180 mm

0,8 mm
1/1
MF-137T
180 mm

0,4 mm
1/1
MF-138T
210 mm

0,4 mm
1/1
MF-139T
210 mm

0,8 mm
1/1
MF-140T
210 mm

0,8 mm
1/1
MF-141T
210 mm

MF-130T/MF-141T

Mikro-Knüpfpinzette
Ø 8 mm



6 x 0,4 mm
1/1
MF-362T
150 mm

6 x 0,4 mm
1/1
MF-363T
150 mm

6 x 0,8 mm
1/1
MF-364T
180 mm

6 x 0,8 mm
1/1
MF-365T
180 mm

6 x 0,8 mm
1/1
MF-366T
210 mm

6 x 0,8 mm
1/1
MF-367T
210 mm


MF-362T/MF-367T

Titanium Ring Forceps

Smooth

Serrated

Mikro-Ring Forceps
Ø 8 mm



1 x 0,5 mm



1/1

MF-142T
150 mm

2 x 1 mm



1/1

MF-145T
150 mm

1 x 0,5 mm



1/1

MF-342T
150 mm

2 x 1 mm



1/1

MF-345T
150 mm

1 x 0,5 mm



1/1

MF-143T
180 mm

2 x 1 mm



1/1

MF-146T
180 mm

1 x 0,5 mm



1/1

MF-343T
180 mm

2 x 1 mm



1/1

MF-346T
180 mm

1 x 0,5 mm



1/1

MF-144T
210 mm

2 x 1 mm



1/1

MF-147T
210 mm

1 x 0,5 mm



1/1

MF-344T
210 mm

2 x 1 mm




1/1

MF-347T
210 mm

MF-142T/MF-347T

Titanium Tissue Forceps

Atraumatic Micro Forcep



Atrauma-Fino
1,2 mm



MF-348T
150 mm

Atrauma-Fino
1,2 mm



MF-349T
150 mm

Atrauma-Standard
1,8 mm



MF-354T
150 mm

Atrauma-Standard
1,8 mm



MF-355T
150 mm

Atrauma-Fino
1,2 mm



MF-350T
180 mm

Atrauma-Fino
1,2 mm



MF-351T
180 mm

Atrauma-Standard
1,8 mm



MF-356T
180 mm

Atrauma-Standard
1,8 mm



MF-357T
180 mm

Atrauma-Fino
1,2 mm



MF-352T
210 mm

Atrauma-Fino
1,2 mm



MF-353T
210 mm

Atrauma-Standard
1,8 mm



MF-358T
210 mm


Atrauma-Standard
1,8 mm




MF-359T
210 mm

MF-348T/MF-359T


Titanium Tissue Forceps




DE BAKEY
BD-80T/BD-163T




BD-160T
160 mm




BD-161T
180 mm




BD-162T
200 mm




BD-163T
240 mm




BD-80T
160 mm




BD-81T
180 mm



BD-84T
200 mm



BD-86T
240 mm



DE BAKEY
BD-82T
160 mm


Micro Scissors & Needle Holders

Approximately \varnothing 8 mm




Micro Needle Holders


Micro Scissors



FD-562T FD-563T



FD-574T FD-575T




FD-587T FD-588T

FD-562T/FD-563T
180 mm


FD-574T/FD-575T
180 mm

FD-587T/FD-588T
180 mm


Fine Artery Forceps




HALSTED-MOSQUITO
P-401/P-402
120 mm



HALSTED-MOSQUITO
P-403/P-404
120 mm



KELLY
P-405/P-406
140 mm



CRILE
P-407/P-408
140 mm

Ceramic Instruments



Castroveijo Caliper
0042
20mm



0043
Castroveijo Scissor
10cm CVD



0047
Trout Barraquer
Needle Holder
13cm CVD



0048
Vannas Scissor
8cm CVD



0062-
Castroviejo Scissor
CVD 10cm



0069
Westcott Tenotomy Scissor
10cm CVD



0049
McPherson Time Forcep
10cm



0051
Harms Colibri Forcep
0.12mm 7cm



0053
MicroBishop Tissue Forcep
8.5cm



0056
Utrata Capsular Hexis Forcep
Moria 11cm



0057
Suture time Forcep
Straight 7mm

Ceramic Instruments



0060
Suture Time forcep
CVD 7mm



0064
Castroviejo Suture Forcep
Straight 0.12mm



0065
Castroviejo Suture Forcep
Straight 0.12mm



0067
Utrata Capsular Hexis Forcep
11cm



0070- IOL Hook



0066
Barraquer Eye Speculum
4cm



INDEX

A

Arrowsmith.....	03
Anis.....	13,16,21,28,32,52
Azar.....	24
Agricola.....	40
Arruga.....	41
Air Injection.....	44
Aldgerbrush.....	45
Ambrosio.....	46
Air Injection Cannula,30ga.....	54
Arcuate Incision Marker.....	59

B

Botvin.....	01
Bross.....	02
Bishop-Harmom.....	05
Barken.....	04
Barraquer.....	07,18,25,31,34,32,33
Biard.....	08
Backhaus.....	10
Baby jone.....	10
Berke.....	11
Bechert.....	19
Bonn.....	26
Bunge.....	24
Bard-Parker Handel.....	48
Bard-Parker Blade.....	48
Bishop-Harmon.....	50
Biploar Forceps.....	51
Biploar.....	51
Biplor Cord.....	51

C

Colibri.....	01,02
Castroviejo.....	12,13,14,20,30,31,36
Culler.....	18
Clayman.....	16,19
Cook.....	35
Castroviejo.....	45
Chuck Handle.....	48
Cannulas.....	52,53
Catrract Set.....	54

D

Dyson-Colibri.....	01
Dressing.....	04
Desmarres.....	08
Drysdale.....	19
Drews.....	31,49,52

E

Elsching.....	06,20
Erhardt.....	11
Eye Scissors.....	17
Emery.....	42
Elliot.....	45
Enucleation Set.....	55,59

F

Farkas-Bracken.....	06
Fine.....	07

Francis.....	08,43
Fixation Forceps.....	28
Fenzl.....	29
Feaster-Barraquer.....	34
Flieringe.....	39
Finskind-castroviejo.....	36

G

Gill.....	01,21,43
Gills-Colibri.....	01
Girards.....	03
Gradle.....	07
Graefe.....	06
Gills-Vannas.....	15
Gimbel.....	18
Gills-Welsh.....	23
Girards.....	22
Gass.....	21
Graether.....	26,24,26,47
Galand.....	28
Grandon-Barraquer.....	34
Guyton-Park.....	38
Golf-Club.....	43
Girard-swan.....	42

H

Harms-Colibri.....	01
Hofmann-polack.....	02
Harms.....	07
Hunt.....	08
Hartman.....	10
Halsted.....	10
Hirschman.....	22,27
Hydro-Dissection.....	49
Hydrodissection.....	50

J

Jewelers.....	09
Jameson.....	11
Jaffe-Bechert.....	19
Jaffe.....	22,30,41
Jaffe-Knolle.....	26
Jaffe-Maltzman.....	28
Jaeger.....	41
J-Shaped Cannula.....	50
Jensen.....	49

K


Kuhnt.....	06
Katena.....	07,05
Katena-Vannas.....	15
Knolle.....	18,26
Kraff.....	16
Kimura.....	20
Kelly.....	21
Kansas.....	23,25
Keener-Arit.....	23
Keener.....	23
Kuglen.....	29,38

INDEX



Kalt.....	31	Radil Keratotomy Set.....	56
Kartz-Barraquer.....	35	S	
Knapp.....	40	Sauer.....	05
Kelman.....	44	Suture Removal Forceps.....	09
Ketena.....	45	Serrefine.....	10
Knolle-Kelman.....	47	Shepard II.....	11
Kartz.....	50,49	Shaepard.....	16
Kraff.....	49	Shepard.....	26
Knolle-Pearce.....	52	Sinsky.....	27
L		Sauer.....	35
Lester.....	04,23	Smirmaul.....	39
Lambert.....	08	Schott.....	39
Lordan.....	08	Stevenson.....	40
Lindstrom.....	13	Stiletto Blade.....	42
Lewis.....	20	Spud & Gouge.....	43
Lieppman.....	22,47	Simcoe.....	46
Lewcky.....	27	Swiss.....	46
Lancaster.....	37	Swiss Model.....	45
Lester-Burcs.....	37	Stiletto Knife.....	48
Lamellar Blade.....	48	Simcoe.....	49
Lid set.....	57	Silicone.....	51
Lacrimal Set.....	58	Sheet.....	52
M		Simcoe.....	52
Maumenee-Colibri.....	03	Sterilizing Case.....	58
Manhattan Eye & Ear.....	05	Sterilizing Case.....	59
McPherson-Vannas.....	15	T	
McPherson-Westcott.....	14	Tennant-Coilibri.....	03
Morrison.....	23	Tissue.....	04
Meyerhoefer.....	21	Tannant.....	07,19
McPherson.....	30,36	Tennant-Troutman.....	06
Maltzman-Fenzl.....	28	Troutman.....	06,33,44
McInyre.....	32	Troutman-Castroviejo.....	12
McIntyer.....	34	Twist-Hook.....	28
Meaumennee-Park.....	38	Tennant.....	35
McIntrya.....	47	Tooke.....	43
Miniature Blade.....	48	Thurmond.....	53
Mcintyre-Binkhorst.....	52	V	
N		Vannas.....	15
New Orleans.....	20	W	
Nightingale.....	21	Wies.....	06
Nishi.....	53	Westcott.....	14
O		Wilder.....	20
Osher.....	13	Williams.....	37
P		Wills Eye Hospital.....	43
Polack.....	02	West.....	44
Pierse.....	02		
Pisacano.....	18		
Paton.....	25		
R			
Rubman Bores.....	03		
Rentsch"Boat Hook".....	29		
Reversed Simcoe.....	46		
Randolph.....	44		
Rubman-Katzin.....	49		

Specialized In Manufacturing of Eye Instruments



QUALITY INSTRUMENTS!

MODREN RANGE OF EYE INSTRUMENTS